


# Easier orientation, navigation and wide-angle survey

## ZEISS Fish Eye Mode



Sample image, using the Fish Eye license.

The Fish Eye license presents an extremely large field of view by permitting fish eye distortions of the image. In this mode the usual measurement tools are disabled. Choosing Fish Eye Mode will select the necessary operating conditions for the column that are required. This will include changing the accelerating voltage.

### Highlights

#### Wide-angle view

Benefit from image samples or entire sample holders with a very large field of view

#### Easier orientation and navigation

Find your region of interest (ROI) even on large samples

### Upgrade your ZEISS microscope

Upgrade your ZEISS microscope with the Fish Eye license and benefit from easier orientation, navigation and a wide-angle view.

Contact us to learn more about the Fish Eye license and how your processes will benefit from an upgrade:  
[microscopy@zeiss.com](mailto:microscopy@zeiss.com)

### Availability

- Crossbeam series
- Sigma series
- EVO series
- AURIGA series
- ULTRA series
- SUPRA series
- NVision series
- NEON series
- 15xx XB series
- 15xx series

Not for therapeutic use, treatment or medical diagnostic evidence. Not all products are available in every country. Contact your local ZEISS representative for more information.  
EN\_43\_012\_068 | CZ-03/2020 | Design, scope of delivery and technical progress subject to change without notice. | © Carl Zeiss Microscopy GmbH



[microscopy@zeiss.com](mailto:microscopy@zeiss.com)  
[www.zeiss.com/microscopy](http://www.zeiss.com/microscopy)



Seeing beyond